

Product datasheet for **TA398346**

Fab Sheep IgG (H&L) Antibody Biotin Conjugated

Product data:

| | |
|------------------------|---|
| Product Type: | Secondary Antibodies |
| Product Name: | Fab Sheep IgG (H&L) Antibody Biotin Conjugated |
| Applications: | ELISA, IHC, WB |
| Recommended Dilution: | WB: 1:2,000 - 1:10,000 IHC: 1:1,000 - 1:5,000 ELISA: 1:20,000 - 1:100,000 |
| Host: | Rabbit |
| Immunogen: | Sheep IgG whole molecule |
| Formulation: | 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2 |
| Reconstitution Method: | Restore with deionized water (or equivalent) - Reconstitution Volume: 500 µL |
| Concentration: | 0.5 mg/mL - lot specific |
| Conjugation: | Biotin |
| Storage: | Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. |
| Note: | Suitable for immunoblotting, ELISA, immunohistochemistry, immunomicroscopy as well as other antibody based assays using streptavidin or avidin conjugates requiring extremely low background levels, absence of F(c) mediated binding, lot-to-lot consistency, high titer and specificity. |



[View online »](#)