

## Product datasheet for **TA398360**

### Ighg3 Rabbit Monoclonal Antibody [Clone ID: RM218]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	RM218
Applications:	ELISA, WB
Recommended Dilution:	ELISA: 0.005 ug/mL-0.2 ug/mL WB (WB): 0.1 ug/mL-0.5 ug/mL
Reactivity:	Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Mouse IgG3
Specificity:	This antibody reacts to the Fab region of mouse IgG3. No cross reactivity with mouse IgG1, IgG2a, IgG2b, IgM, IgA, IgE, human IgG, or rat IgG.
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	1 mg/mL - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	Immunoglobulin heavy constant gamma 3
Database Link:	<a href="#">Entrez Gene 380795 Mouse P03987</a>



[View online »](#)