

Product datasheet for **TA398370**

H3FA (HIST1H3A) Rabbit Monoclonal Antibody [Clone ID: RM157]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM157
Applications:	ELISA, LMNX, WB
Recommended Dilution:	WB: 0.2 ug/mL – 1 ug/mL ELISA: 0.1 ug/mL - 0.5 ug/mL Luminex: 0.1 ug/mL – 0.5 ug/mL
Reactivity:	Human, Vertebrates
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A trimethyl-peptide corresponding to Trimethyl-Histone H3 (Lys79)
Specificity:	This antibody reacts to Histone H3 trimethylated at Lysine 79 (K79me3). No cross reactivity with monomethylated Lysine 79 (K79me1) or dimethylated Lysine 79 (K79me2), or other methylations in histone H3.
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	1 mg/mL - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	histone cluster 1, H3a
Database Link:	Entrez Gene 8350 Human P84243
Synonyms:	H3/A; H3/b; H3/c; H3/d; H3/f; H3/h; H3/i; H3/j; H3/k; H3/l; H3FA; H3FB; H3FC; H3FD; H3FE; H3FH; H3FI; H3FJ; H3FK; H3FL



[View online »](#)