

## Product datasheet for **TA398374**

### **H3FA (HIST1H3A) Rabbit Monoclonal Antibody [Clone ID: RM322]**

#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	RM322
Applications:	LMNX, WB
Recommended Dilution:	WB: 0.1 ug/mL – 1 ug/mL Luminex: 0.01 ug/mL – 1 ug/mL
Reactivity:	Human, Vertebrates
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A dimethyl-peptide corresponding to dimethyl-Histone H3 (Lys9) and an acetyl-peptide corresponding to Acetyl-Histone H3 (Lys14)
Specificity:	This antibody reacts to human Dimethyl-Histone H3 (Lys9) and Acetyl-Histone H3 (Lys14)
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	1 mg/mL - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	histone cluster 1, H3a
Database Link:	<a href="#">Entrez Gene 8350 Human P84243</a>
Synonyms:	H3/A; H3/b; H3/c; H3/d; H3/f; H3/h; H3/i; H3/j; H3/k; H3/l; H3FA; H3FB; H3FC; H3FD; H3FF; H3FH; H3FI; H3FJ; H3FK; H3FL



[View online »](#)