

Product datasheet for **TA398382**

H3FA (HIST1H3A) Rabbit Monoclonal Antibody [Clone ID: RM186]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM186
Applications:	ELISA, ICC, LMNX, WB
Recommended Dilution:	Western Blot: 1 ug/mL – 2 ug/mL ELISA: 0.5 ug/mL - 1 ug/mL Luminex: 0.5 ug/mL – 2 ug/mL Immunocytochemistry: 0.5 ug/mL- 2 ug/mL
Reactivity:	Human, Vertebrates
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	An unmodified peptide corresponding to the N-terminus of Histone H3.
Specificity:	This antibody recognizes Histone H3 that is unmodified at Lys4 but does not recognize acetylated or methylated Lys4. The antibody binding specificity allows for modifications of Arg2, Thr3, and/or other modifications in Histone H3.
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	1 mg/mL - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	histone cluster 1, H3a
Database Link:	Entrez Gene 8350 Human Q16695
Synonyms:	H3/A; H3/b; H3/c; H3/d; H3/f; H3/h; H3/i; H3/j; H3/k; H3/l; H3FA; H3FB; H3FC; H3FD; H3FF; H3FH; H3FI; H3FJ; H3FK; H3FL



[View online »](#)