

Product datasheet for **TA398393**

H2A.Z (H2AFZ) Rabbit Monoclonal Antibody [Clone ID: RM221]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM221
Applications:	ELISA, ICC, LMNX, WB
Recommended Dilution:	WB: 0.5 ug/mL – 2 ug/mL ELISA: 0.2 ug/mL - 1 ug/mL Luminex: 0.05 ug/mL – 0.5 ug/ml Immunocytochemistry: 1 ug/mL- 2 ug/mL
Reactivity:	Human, Vertebrates
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	An acetyl-peptide corresponding to Acetyl-Histone H2A.Z (Lys4)
Specificity:	This antibody reacts to Histone H2A.Z acetylated at Lysine 4 (K4ac). No cross reactivity with non-modified Lysine 4 or other acetylated Lysines in histone H2A
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	1 mg/mL - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	H2A histone family member Z
Database Link:	Entrez Gene 3015 Human P0C0S5
Synonyms:	H2A.z; H2A/z; H2AZ; MGC117173



[View online »](#)