

Product datasheet for **TA398404**

H2BC3 Rabbit Monoclonal Antibody [Clone ID: RM455]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM455
Applications:	ELISA, LMNX, WB
Recommended Dilution:	WB: 0.01 ug/mL - 1 ug/mL ELISA: 0.2 ug/mL - 1 ug/mL Luminex: 0.01 ug/mL - 0.5 ug/mL
Reactivity:	Human, Vertebrates
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	An acetyl-peptide corresponding to Acetyl-Histone H2B (Lys5).
Specificity:	This antibody reacts to Histone H2B acetylated at Lysine 5 (K5ac). No cross reactivity with other acetylated Lysines in histones.
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	1 mg/mL - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	histone cluster 1, H2bb
Database Link:	Entrez Gene 3018 Human P33778
Synonyms:	H2B.1; H2B.f; H2B/f; H2BFF; MGC119804



[View online »](#)