

## Product datasheet for **TA398412**

### Rabbit Monoclonal Antibody [Clone ID: RMH04]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	RMH04
Recommended Dilution:	ELISA: 50ng/well – 200ng/wel
Reactivity:	Ape
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	25-OH Vitamin D3 BSA conjugate
Specificity:	This antibody reacts to 25-OH Vitamin D3. It does not cross-react with Vitamin D3, Vitamin D2, and 25-OH Vitamin D2.
Formulation:	50% Glycerol/PBS and 0.09% sodium azide
Concentration:	1 mg/mL - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.



[View online »](#)