

Product datasheet for **TA398413**

Rabbit Monoclonal Antibody [Clone ID: RM343]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM343
Recommended Dilution:	Competitive ELISA: 5 ng/well – 10 ng/well on a 96-well plate pre-coated with a goat anti-rabbit IgG polyclonal antibody
Reactivity:	Ape
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Estradiol BSA Conjugate.
Specificity:	This antibody specifically reacts to Estradiol. It does not cross-react with Estradiol 17 β , Estrone, Estriol, Progesterone, or DHEA-SO ₄ .
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.



[View online »](#)