

## Product datasheet for **TA398414**

### Rabbit Monoclonal Antibody [Clone ID: RM434]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	RM434
Recommended Dilution:	Competitive ELISA: 1 ng/well – 10 ng/wel
Reactivity:	Ape
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Progesterone-KLH
Specificity:	This antibody reacts to Progesterone
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.



[View online »](#)