

## Product datasheet for **TA398416**

### Rabbit Monoclonal Antibody [Clone ID: RM466]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	RM466
Applications:	ELISA
Recommended Dilution:	ELISA: 0.05 ug/mL – 0.5 ug/m
Reactivity:	Ape
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	KLH-conjugated cAMP
Specificity:	This antibody reacts to cAMP (Cyclic Adenosine Monophosphate). No cross reactivity with other cyclic nucleotide or nucleoside phosphate.
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	1 mg/mL - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.



[View online »](#)