

## Product datasheet for **TA398422**

### Rabbit Monoclonal Antibody [Clone ID: RM477]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	RM477
Applications:	Dot, ELISA
Recommended Dilution:	Dot Blot: 0.5 ug/mL - 2 ug/mL ELISA: 0.5 ug/mL - 2 ug/m
Reactivity:	Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	KLH-conjugated 5-Formylcytosines
Specificity:	This antibody reacts to 5-Formylcytosines (5-fC) in DNA. No cross reactivity with non-modified cytosine and other modified cytosines in DNA.
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	1 mg/mL - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.



[View online »](#)