

Product datasheet for **TA398434**

RSK1 p90 (RPS6KA1) Rabbit Monoclonal Antibody [Clone ID: RM233]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM233
Applications:	IHC, WB
Recommended Dilution:	Immunohistochemistry (IHC): 1:200-1:500 dilution WB (WB): 1:1000-1:2000 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A phospho-peptide corresponding to human phospho-Rsk1 (Thr359/Ser363)
Specificity:	This antibody reacts to human RSK1 only when phosphorylated at Thr359/Ser363.
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	ribosomal protein S6 kinase A1
Database Link:	Entrez Gene 6195 Human Q15418
Synonyms:	HU-1; MAPKAPK1A; p90S6K; pp90RSK1; RSK; RSK-1; RSK1



[View online »](#)