

Product datasheet for **TA398442**

EGFR Rabbit Monoclonal Antibody [Clone ID: RM443]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM443
Applications:	WB
Recommended Dilution:	WB:1:1000-1:400
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A phospho-peptide corresponding to human Phospho- EGFR (Y1068)
Specificity:	This antibody reacts to human EGF Receptor only when phosphorylated at Y1068
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	epidermal growth factor receptor
Database Link:	Entrez Gene 1956 Human P00533
Synonyms:	ERBB; ERBB1; HER1; mENA; PIG61



[View online »](#)