

Product datasheet for **TA398446**

Rabbit Monoclonal Antibody [Clone ID: RM454]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM454
Applications:	WB
Recommended Dilution:	WB: 1:1000 -1:3000 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A phospho-peptide corresponding to Phospho-Aurora A (Thr288)/Aurora B (Thr232)/Aurora C (Thr198)
Specificity:	This antibody reacts to human Aurora A/B/C only when phosphorylated at either Thr288, Thr232 or Thr198 respectively. There is no cross-reactivity to Aurora A/B/C that is not phosphorylated.
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Database Link:	Q9UQB9 , Q96GD4 , O14965

[View online »](#)