

Product datasheet for **TA398447**

Tau (MAPT) Rabbit Monoclonal Antibody [Clone ID: RM461]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM461
Applications:	IHC, WB
Recommended Dilution:	IHC: 1:50 – 1:200 dilution WB: 1:1000 -1:2000 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A phospho-peptide corresponding to human Tau (Ser396)
Specificity:	This antibody reacts to Tau only when phosphorylated at Ser396.
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	microtubule associated protein tau
Database Link:	Entrez Gene 4137 Human P10636
Synonyms:	DDPAC; FLJ31424; FTDP-17; MAPTL; MGC138549; MSTD; MTBT1; MTBT2; PHF-tau; PPND; TAU



[View online »](#)