

Product datasheet for **TA398448**

Rabbit Monoclonal Antibody [Clone ID: RM464]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM464
Applications:	IHC, WB
Recommended Dilution:	IHC: 1:100 - 1:200 dilution WB: 1:200 -1:2000 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A phospho-peptide corresponding to human Phospho- JNK1/2/3 (Thr183/Tyr185)
Specificity:	This antibody reacts to human JNK1, JNK2, and JNK3 only when phosphorylated at Thr183 and Tyr185.
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Database Link:	<u>P45983/P45984/P53779</u>

[View online »](#)