

## Product datasheet for **TA398449**

### Tubulin (TUBA1B) Rabbit Monoclonal Antibody [Clone ID: RM318]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	RM318
Applications:	IHC, WB
Recommended Dilution:	Immunocytochemistry (ICC): 1:1000-1:5000 dilution WB (WB): 1:1000-1:2000 dilution
Reactivity:	Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	An acetyl-peptide corresponding to Acetyl- alpha-Tubulin (Lys40)
Specificity:	This antibody reacts to alpha-Tubulin acetylated at Lysine 40. No cross reactivity to non-acetylated alpha-Tubulin at Lysine 40. This antibody may also react to mouse or rat Acetyl-alpha-Tubulin at Lysine 40, as predicted by immunogen homology.
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	tubulin alpha 1b
Database Link:	<a href="#">P68363</a>
Synonyms:	K-ALPHA-1



[View online »](#)