

Product datasheet for **TA398450**

Cytokeratin 5 (KRT5) Rabbit Monoclonal Antibody [Clone ID: RM226]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM226
Applications:	IHC, WB
Recommended Dilution:	Immunohistochemistry (IHC): 1:100 dilution WB (WB): 1:1000 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to the C-terminus of Human Cytokeratin 5
Specificity:	This antibody reacts to Human Cytokeratin 5. This antibody may also react to bovine, mouse, or rat Cytokeratin, as predicted by immunogen homology.
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	keratin 5
Database Link:	Entrez Gene 3852 Human P13647
Synonyms:	CK-5; CK5; Cytokeratin-5; DDD; EBS2; K5; Keratin-5; KRT5A



[View online »](#)