

Product datasheet for **TA398451**

Her2 (ERBB2) Rabbit Monoclonal Antibody [Clone ID: RM228]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM228
Applications:	IHC, WB
Recommended Dilution:	Immunohistochemistry (IHC): 1:100 dilution WB (WB): 1:1000 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to the C-terminus of Human c-erbB-2/HER-2
Specificity:	This antibody reacts to Human c-erbB-2/HER-2. This antibody may also react to mouse or rat c-erbB-2/HER-2, as predicted by immunogen homology.
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	erb-b2 receptor tyrosine kinase 2
Database Link:	Entrez Gene 2064 Human P04626
Synonyms:	C-erbB-2; CD340; erbB-2; HER-2; HER-2/neu; HER2; herstatin; NEU; NGL; p185erbB2; TKR1



[View online »](#)