

Product datasheet for **TA398453**

Desmin (DES) Rabbit Monoclonal Antibody [Clone ID: RM234]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM234
Applications:	IHC, WB
Recommended Dilution:	Immunohistochemistry (IHC): 1:1000-1:4000 dilution WB (WB): 1:1000-1:2000 dilution
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to human Desmin
Specificity:	This antibody reacts to human and mouse Desmin. This antibody may also react to rat and Bovine Desmin, as predicted by immunogen homology.
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	desmin
Database Link:	Entrez Gene 1674 Human P17661
Synonyms:	CMD1I; CSM1; CSM2; desmin; FLJ12025; FLJ39719; FLJ41013; FLJ41793



[View online »](#)