

Product datasheet for **TA398454**

Cyclin D1 (CCND1) Rabbit Monoclonal Antibody [Clone ID: RM241]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM241
Applications:	IHC, WB
Recommended Dilution:	Immunohistochemistry (IHC): 1:500-1:1000 dilution WB (WB): 1:1000-1:2000 dilution
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to Cyclin D1
Specificity:	This antibody reacts to human Cyclin D1. This antibody may also react to bovine, mouse or rat Cyclin D1, as predicted by immunogen homology.
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	cyclin D1
Database Link:	Entrez Gene 595 Human P24385
Synonyms:	BCL1; D11S287E; PRAD1; U21B31



[View online »](#)