

## Product datasheet for **TA398458**

### GFAP Rabbit Monoclonal Antibody [Clone ID: RM246]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	RM246
Applications:	IHC, WB
Recommended Dilution:	Immunohistochemistry (IHC): 1:200-1:500 dilution WB (WB): 1:1000-1:3000 dilution
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to the N-terminus of human GFAP
Specificity:	This antibody reacts to human and mouse GFAP (Glial fibrillary acidic protein).
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	glial fibrillary acidic protein
Database Link:	<a href="#">Entrez Gene 2670 Human P14136</a>
Synonyms:	FLJ45472



[View online »](#)