

Product datasheet for **TA398466**

Calponin (CNN1) Rabbit Monoclonal Antibody [Clone ID: RM262]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM262
Applications:	ICC, WB
Recommended Dilution:	WB (WB): 1:1000-1:2000 dilution Immunocytochemistry (ICC): 1:100-1:200
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to the C-terminus of human Calponin-1
Specificity:	This antibody reacts to human or mouse Calponin-1
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	calponin 1
Database Link:	Entrez Gene 1264 Human P51911
Synonyms:	Sm-Calp; SMCC



[View online »](#)