

Product datasheet for **TA398467**

N Cadherin (CDH2) Rabbit Monoclonal Antibody [Clone ID: RM259]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM259
Applications:	WB
Recommended Dilution:	WB (WB): 1:1000-1:2000 dilutio
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to Human N-Cadherin
Specificity:	This antibody reacts to human N-cadherin. This antibody may also react to mouse or rat N-cadherin, as predicted by immunogen homology.
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	cadherin 2
Database Link:	Entrez Gene 1000 Human P19022
Synonyms:	CD325; CDHN; CDw325; N-cadherin; NCAD



[View online »](#)