

Product datasheet for **TA398468**

Paxillin (PXN) Rabbit Monoclonal Antibody [Clone ID: RM256]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM256
Applications:	ICC, IHC, WB
Recommended Dilution:	Immunocytochemistry (ICC): 1:1000-1:2000 dilution WB (WB): 1:1000-1:2000 dilution Immunohistochemistry (IHC): 1:1000 dilution
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to the N-terminus of human Paxillin
Specificity:	This antibody reacts to human Paxillin. This antibody may also react to mouse or rat Paxillin, as predicted by immunogen homology.
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	paxillin
Database Link:	Entrez Gene 5829 Human P49023
Synonyms:	FLJ16691; paxillin



[View online »](#)