

Product datasheet for **TA398471**

CD44 Rabbit Monoclonal Antibody [Clone ID: RM264]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM264
Applications:	IHC, WB
Recommended Dilution:	Immunohistochemistry (IHC): 1:1000-1:5000 dilution WB (WB): 1:1000-1:2000 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to the extracellular domain of human CD44.
Specificity:	This antibody reacts to the extracellular domain of human CD44.
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	CD44 molecule (Indian blood group)
Database Link:	Entrez Gene 960 Human P16070
Synonyms:	CD44R; CDW44; CSPG8; ECMR-III; Epican; HCELL; HUTCH-I; IN; LHR; MC56; MDU2; MDU3; MGC10468; MIC4; MUTCH-I; PGP-1; PGP-I; Pgp1



[View online »](#)