

Product datasheet for **TA398475**

Integrin alpha 4 (ITGA4) Rabbit Monoclonal Antibody [Clone ID: RM268]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM268
Applications:	IHC, WB
Recommended Dilution:	Immunohistochemistry (IHC): 1:1000-1:5000 dilution WB (WB): 1:1000-1:2000 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to the C-terminus of human Integrin alpha 4 (ITGA4).
Specificity:	This antibody reacts to the cytoplasmic domain of human Integrin alpha 4 (ITGA4).
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	integrin subunit alpha 4
Database Link:	Entrez Gene 3676 Human P13612
Synonyms:	CD49D; IA4; MGC90518; VLA-4



[View online »](#)