

Product datasheet for **TA398476**

DIABLO Rabbit Monoclonal Antibody [Clone ID: RM271]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM271
Applications:	IHC, WB
Recommended Dilution:	Immunohistochemistry (IHC): 1:500-1:1000 dilution WB (WB): 1:1000-1:2000 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to the C-terminus of human Smac (Diablo homolog, mitochondrial).
Specificity:	This antibody reacts to human Smac (Diablo homolog, mitochondrial).
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	diablo IAP-binding mitochondrial protein
Database Link:	Entrez Gene 56616 Human Q9NR28
Synonyms:	0610041G12Rik; DIABLO-S; FLJ10537; FLJ25049; SMAC; SMAC3



[View online »](#)