

Product datasheet for **TA398483**

Aurora B (AURKB) Rabbit Monoclonal Antibody [Clone ID: RM278]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM278
Applications:	IHC, WB
Recommended Dilution:	Immunohistochemistry (IHC): 1:200-1:500 dilution WB (WB): 1:1000-1:2000 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to the N-terminus of human Aurora kinase B.
Specificity:	This antibody reacts to human Aurora kinase B.
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	aurora kinase B
Database Link:	Entrez Gene 9212 Human Q96GD4
Synonyms:	AIK2; AIM-1; AIM1; ARK2; AurB; aurkb-sv1; aurkb-sv2; aurora-1; aurora-B; IPL1; STK-1; STK5; STK12



[View online »](#)