

## Product datasheet for **TA398490**

### EGFR Rabbit Monoclonal Antibody [Clone ID: RM294]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	RM294
Applications:	IHC, WB
Recommended Dilution:	Immunohistochemistry (IHC): 1:100-1:200 dilution WB (WB): 1:200-1:500 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to the C-terminus of human EGF Receptor (c-ErbB-1).
Specificity:	This antibody reacts to human EGF Receptor (c-ErbB-1). No cross reactivity with HER2, HER3, or HER4
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	epidermal growth factor receptor
Database Link:	<a href="#">Entrez Gene 1956 Human P00533</a>
Synonyms:	ERBB; ERBB1; HER1; mENA; PIG61



[View online »](#)