

Product datasheet for **TA398494**

alpha Actinin 4 (ACTN4) Rabbit Monoclonal Antibody [Clone ID: RM287]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM287
Applications:	IHC, WB
Recommended Dilution:	Immunohistochemistry (IHC): 1:100-1:200 dilution WB (WB): 1:1000-1:2000 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to the N-terminus of human alpha-Actinin-4
Specificity:	This antibody reacts to human alpha-Actinin-4.
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	actinin alpha 4
Database Link:	Entrez Gene 81 Human O43707
Synonyms:	ACTININ-4; DKFZp686K23158; FSGS; FSGS1



[View online »](#)