

Product datasheet for **TA398495**

PSD95 (DLG4) Rabbit Monoclonal Antibody [Clone ID: RM288]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM288
Applications:	IHC, WB
Recommended Dilution:	Immunohistochemistry (IHC): 1:100-1:200 dilution WB (WB): 1:1000-1:8000 dilution
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to the N-terminus of human PSD95
Specificity:	This antibody reacts to human and mouse PSD95. It may also react to rat PSD95, as predicted by immunogen homology
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	discs large MAGUK scaffold protein 4
Database Link:	Entrez Gene 1742 Human P78352
Synonyms:	FLJ97752; FLJ98574; PSD-95; PSD95; SAP-90; SAP90



[View online »](#)