

## Product datasheet for **TA398497**

### CD11b (ITGAM) Rabbit Monoclonal Antibody [Clone ID: RM290]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	RM290
Applications:	IHC, WB
Recommended Dilution:	Immunohistochemistry (IHC): 1:100-1:200 dilution WB (WB): 1:200-1:500 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to the C-terminus of human CD11b (Integrin alpha-M)
Specificity:	This antibody reacts to human CD11b (Integrin alpha-M)
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	integrin subunit alpha M
Database Link:	<a href="#">Entrez Gene 3684 Human P11215</a>
Synonyms:	CD11B; CR3A; MAC-1; MAC1A; MGC117044; MO1A; SLEB6



[View online »](#)