

## Product datasheet for **TA398501**

### AMPK alpha 1 (PRKAA1) Rabbit Monoclonal Antibody [Clone ID: RM301]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	RM301
Applications:	IHC, WB
Recommended Dilution:	Immunohistochemistry (IHC): 1:100-1:200 dilution WB (WB): 1:100-1:200 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to the residues near the C-terminus of human AMPK alpha-1.
Specificity:	This antibody reacts to human AMPK alpha-1 (AMPKA1) (5'-AMP-activated protein kinase catalytic subunit alpha-1).
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	protein kinase AMP-activated catalytic subunit alpha 1
Database Link:	<a href="#">Entrez Gene 5562 Human Q13131</a>
Synonyms:	AMPK; AMPK1; AMPKa1; MGC33776; MGC57364



[View online »](#)