

Product datasheet for **TA398507**

PD-L1 (Cd274) Rabbit Monoclonal Antibody [Clone ID: RM320]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM320
Applications:	ICC, IHC, WB
Recommended Dilution:	Immunohistochemistry (IHC): 1:200-1:1000 dilution WB (WB): 1:500-1:1000 dilution Immunocytochemistry (ICC): 1:200-1:500 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to the C-terminus of human PD-L1
Specificity:	This antibody reacts to human PD-L1 (CD274) (B7-H1)
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	CD274 antigen
Database Link:	Q9EP73
Synonyms:	B7-H; B7-H1; B7H1; MGC142294; MGC142296; PD-L1; PDCD1L1; PDCD1LG1; PDL1



[View online »](#)