

Product datasheet for **TA398515**

BAG1 Rabbit Monoclonal Antibody [Clone ID: RM310]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM310
Applications:	ICC, IHC, WB
Recommended Dilution:	Immunohistochemistry (IHC): 1:500-1:2000 dilution WB (WB): 1:1000-1:2000 dilution Immunocytochemistry (ICC): 1:200-1:1000 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to residues near the N-terminus of human BAG-1L
Specificity:	This antibody reacts to human BAG-1L. It does not cross react to BAG-1M or BAG-1S.
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	BCL2 associated athanogene 1
Database Link:	Entrez Gene 573 Human Q99933
Synonyms:	BAG-1; HAP; RAP46



[View online »](#)