

## Product datasheet for **TA398516**

### Calretinin (CALB2) Rabbit Monoclonal Antibody [Clone ID: RM324]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	RM324
Applications:	IHC, WB
Recommended Dilution:	Immunohistochemistry (IHC): 1:500-1:1000 dilution WB (WB): 1:2000-1:4000 dilution
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to N-terminus of human Calretinin
Specificity:	This antibody reacts to human Calretinin. It may also react to mouse or rat Calretinin, as predicted by immunogen homology.
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	calbindin 2
Database Link:	<a href="#">Entrez Gene 794 Human P22676</a>
Synonyms:	CAB29; CAL2; calretinin; CR



[View online »](#)