

Product datasheet for **TA398519**

PSMA (FOLH1) Rabbit Monoclonal Antibody [Clone ID: RM327]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM327
Applications:	IHC, WB
Recommended Dilution:	Immunohistochemistry (IHC): 1:500-1:1000 dilution WB (WB): 1:1000-1:2000 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to the C-terminus of human PSMA (Prostate-specific membrane antigen).
Specificity:	This antibody reacts to human PSMA (Prostate-specific membrane antigen).
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	folate hydrolase (prostate-specific membrane antigen) 1
Database Link:	Entrez Gene 2346 Human Q04609
Synonyms:	FGCP; FOLH; GCP2; GCPII; mGCP; NAALAD1; NAALAdase; PSM; PSMA



[View online »](#)