

Product datasheet for **TA398522**

PAX5 Rabbit Monoclonal Antibody [Clone ID: RM331]

Product data:

| | |
|-----------------------|--|
| Product Type: | Primary Antibodies |
| Clone Name: | RM331 |
| Applications: | IHC, WB |
| Recommended Dilution: | Immunohistochemistry (IHC): 1:500-1:1000 dilution WB (WB): 1:1000-1:2000 dilution |
| Reactivity: | Human, Mouse |
| Host: | Rabbit |
| Isotype: | IgG |
| Clonality: | Monoclonal |
| Immunogen: | A peptide corresponding to C-terminus of human Paired box protein Pax-5 |
| Specificity: | This antibody reacts to human Paired box protein Pax-5. It may also react to mouse Paired box protein Pax-5, as predicted by immunogen homology. |
| Formulation: | 50% Glycerol/PBS with 1% BSA and 0.09% sodium azide |
| Concentration: | N/A - lot specific |
| Purification: | Protein A affinity purified from an animal origin-free culture supernatant |
| Conjugation: | Unconjugated |
| Storage: | Store at -20°C as received. |
| Stability: | Stable for 12 months from date of receipt. |
| Gene Name: | paired box 5 |
| Database Link: | Entrez Gene 5079 Human Q02548 |
| Synonyms: | BSAP |



[View online »](#)