

Product datasheet for **TA398524**

MelanA (MLANA) Rabbit Monoclonal Antibody [Clone ID: RM333]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM333
Applications:	IHC, WB
Recommended Dilution:	Immunohistochemistry (IHC): 1:500-1:1000 dilution WB (WB): 1:500-1:1000 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to the N-terminus of human MART1 (Melanoma antigen recognized by T-cells-1).
Specificity:	This antibody reacts to human MART1 (Melanoma antigen recognized by T-cells-1).
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	melan-A
Database Link:	Entrez Gene 2315 Human Q16655
Synonyms:	MART-1; MART1; melan-A; OTTHUMP00000021037; OTTHUMP00000021038



[View online »](#)