

Product datasheet for **TA398525**

SFTPA1B Rabbit Monoclonal Antibody [Clone ID: RM334]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM334
Applications:	IHC, WB
Recommended Dilution:	Immunohistochemistry (IHC): 1:500-1:1000 dilution WB (WB): 1:1000-1:2000 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to the C-terminus of human Surfactant protein A (SP-A).
Specificity:	This antibody reacts to human Surfactant protein A (SP-A).
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	surfactant, pulmonary-associated protein A1
Database Link:	Q8IWL1



[View online »](#)