

Product datasheet for **TA398534**

GST3 (GSTP1) Rabbit Monoclonal Antibody [Clone ID: RM347]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM347
Applications:	IHC, WB
Recommended Dilution:	Immunohistochemistry (IHC): 1:1000-1:20000 dilution WB (WB): 1:1000-1:2000 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to the C-terminus of human GST3/GST-pi.
Specificity:	This antibody reacts to human GST3/GST-pi (Glutathione S-transferase P).
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	glutathione S-transferase pi 1
Database Link:	Entrez Gene 2950 Human P09211
Synonyms:	DFN7; FAES3; GST3; GSTP; GSTP1-1; PI



[View online »](#)