

Product datasheet for **TA398536**

Stathmin 1 (STMN1) Rabbit Monoclonal Antibody [Clone ID: RM350]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM350
Applications:	IHC, WB
Recommended Dilution:	Immunohistochemistry (IHC): 1:500-1:1000 dilution WB (WB): 1:2000-1:5000 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to the C-terminus of human Stathmin.
Specificity:	This antibody reacts to human Stathmin. It may also react to mouse or rat Stathmin, as predicted by immunogen homology.
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	stathmin 1
Database Link:	Entrez Gene 3925 Human P16949
Synonyms:	C1orf215; FLJ32206; Lag; LAP18; metablastin; MGC138869; MGC138870; OP18; OTTHUMP00000008530; PP17; PP19; PR22; prosolin; SMN



[View online »](#)