

Product datasheet for **TA398540**

Cytokeratin 4 (KRT4) Rabbit Monoclonal Antibody [Clone ID: RM355]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM355
Applications:	IHC, WB
Recommended Dilution:	Immunohistochemistry (IHC): 1:500-1:1000 dilution WB (WB): 1:100-1:1000 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to the C-terminus of human CK-4 (Cytokeratin-4)
Specificity:	This antibody reacts human CK-4 (Cytokeratin-4)
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	keratin 4
Database Link:	Entrez Gene 3851 Human P19013
Synonyms:	CK-4; CK4; CYK4; Cytokeratin-4; FLJ31692; K4; Keratin-4



[View online »](#)