

Product datasheet for **TA398541**

CD117 Rabbit Monoclonal Antibody [Clone ID: RM359]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM359
Applications:	IHC, WB
Recommended Dilution:	Immunohistochemistry (IHC): 1:1000-1:2000 dilution WB (WB): 1:200-1:1000 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to the C-terminus of human CD117/c-Kit
Specificity:	This antibody reacts human CD117/c-Kit
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	KIT proto-oncogene receptor tyrosine kinase
Database Link:	Entrez Gene 3815 Human P10721
Synonyms:	C-Kit; CD117; PBT; SCFR



[View online »](#)