

## Product datasheet for **TA398543**

### Cytokeratin 19 (KRT19) Rabbit Monoclonal Antibody [Clone ID: RM364]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	RM364
Applications:	IHC, WB
Recommended Dilution:	Immunohistochemistry (IHC): 1:500-1:1000 dilution WB (WB): 1:1000-1:2000 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to the C-terminus of human CK-19 (Cytokeratin-19)
Specificity:	This antibody reacts human CK-19 (Cytokeratin-19)
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	keratin 19
Database Link:	<a href="#">Entrez Gene 3880 Human P08727</a>
Synonyms:	CK-19; CK19; Cytokeratin-19; K1CS; K19; Keratin-19; MGC15366



[View online »](#)