

Product datasheet for **TA398544**

SATB2 Rabbit Monoclonal Antibody [Clone ID: RM365]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM365
Applications:	IHC, WB
Recommended Dilution:	Immunohistochemistry (IHC): 1:200-1:1000 dilution WB (WB): 1:1000-1:2000 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to the C-terminus of human DNA-binding protein SATB2
Specificity:	This antibody reacts human DNA-binding protein SATB2
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	SATB homeobox 2
Database Link:	Entrez Gene 23314 Human Q9UPW6
Synonyms:	FLJ21474; FLJ32076; KIAA1034; MGC119474; MGC119477; OTTHUMP00000163655; OTTHUMP00000205774; OTTHUMP00000205775



[View online »](#)