

Product datasheet for **TA398545**

NAPSIN A (NAPSA) Rabbit Monoclonal Antibody [Clone ID: RM366]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM366
Applications:	IHC, WB
Recommended Dilution:	Immunohistochemistry (IHC): 1:500-1:1000 dilution WB (WB): 1:1000-1:2000 dilutio
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to the N-terminus of human Napsin-A
Specificity:	This antibody reacts human Napsin-A
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	napsin A aspartic peptidase
Database Link:	Entrez Gene 9476 Human O96009
Synonyms:	ASP4; KAP; Kdap; NAP1; NAPA; Napsin-1; SNAPA; TA01/TA02



[View online »](#)