

## Product datasheet for **TA398547**

### Thyroid Peroxidase (TPO) Rabbit Monoclonal Antibody [Clone ID: RM368]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	RM368
Applications:	IHC, WB
Recommended Dilution:	Immunohistochemistry (IHC): 1:500-1:1000 dilution WB (WB): 1:200-1:1000 dilution
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to the C-terminus of human TPO (Thyroid peroxidase)
Specificity:	This antibody reacts to human TPO (Thyroid peroxidase)
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	thyroid peroxidase
Database Link:	<a href="#">Entrez Gene 7173 Human P07202</a>
Synonyms:	MSA; OTTHUMP00000115536; TDH2A; thyroperoxidase; TPX



[View online »](#)